



# AC Medical Family Pricelist



# Auguste Cryogenics


Pos	Model	Part No.	Description	Price
1	AC CNL2-S(1)	H-AAA1100128	AC CNL2-S(1), 120 mm canister height 19 mm diameter	14,00 €
2	AC CNL2-S(2)	H-AAA1100129	AC CNL2-S(2), 120 mm canister height 25 mm diameter	15,00 €
3	AC CN L3-S	H-AAA1100130	AC CN L3-S, 120 mm canister height 38 mm diameter	16,00 €
4	AC CN DS3-S	H-AAA1100131	AC CN DS3-S, 120 mm canister height 38 mm diameter	16,00 €
5	AC CN L6-S	H-AAA1100132	AC CN L6-S, 120 mm canister height 38 mm diameter	16,00 €
6	AC CN L10-S	H-AAA1100133	AC CN L10-S, 120 mm canister height 38 mm diameter, used for YDS-10	16,00 €
7	AC CN L15-S	H-AAA1100134	AC CN L15-S, 120 mm canister height 38 mm diameter	16,00 €
8	AC CN L20-S	H-AAA1100135	AC CN L20-S, 120 mm canister height 38 mm diameter, used for YDS-20	16,00 €
9	AC CN L25-S	H-AAA1100136	AC CN L25-S, 120 mm canister height 38 mm diameter	16,00 €
10	AC CN L30-S	H-AAA1100137	AC CN L30-S, 120 mm canister height 38 mm diameter	16,00 €
11	AC CN L35-S	H-AAA1100138	AC CN L35-S, 120 mm canister height 38 mm diameter	16,00 €
12	AC CN L50-S	H-AAA1100139	AC CN L50-S, 120 mm canister height 38 mm diameter	16,00 €
13	AC CN DS6-S	H-AAA1100140	AC CN DS6-S, 120 mm canister height 63 mm diameter	18,00 €
14	AC CN DS8-S	H-AAA1100141	AC CN DS8-S, 120 mm canister height 63 mm diameter	18,00 €
15	AC CN XL10-S	H-AAA1100142	AC CN XL10-S, 120 mm canister height 63 mm diameter	18,00 €
16	AC CN XL20-S	H-AAA1100143	AC CN XL20-S, 120 mm canister height 63 mm diameter	18,00 €
17	AC CN XL30-S	H-AAA1100144	AC CN XL30-S, 120 mm canister height 63 mm diameter	18,00 €
18	AC CN XL35-S	H-AAA1100145	AC CN XL35-S, 120 mm canister height 63 mm diameter	18,00 €
19	AC CN XL50-S	H-AAA1100146	AC CN XL50-S, 120 mm canister height 63 mm diameter	18,00 €
20	AC LL L(S1)	H-AAB1100176	AC LL L(S1), black lockable lid - L (30 mm) - 2 L	19,00 €
21	AC LL L(S2)	H-AAB1100177	AC LL L(S2), black lockable lid - L (35 mm) - 2 L	19,00 €
22	AC CN L13-L	H-AAA1100150	AC CN L13-L, 276 mm canister height 38 mm diameter	20,00 €
23	AC CN L20-L	H-AAA1100151	AC CN L20-L, 276 mm canister height 38 mm diameter	20,00 €
24	AC CN L25-L	H-AAA1100152	AC CN L25-L, 276 mm canister height 38 mm diameter	20,00 €
25	AC CN L30-L	H-AAA1100153	AC CN L30-L, 276 mm canister height 38 mm diameter	20,00 €
26	AC CN L35-L	H-AAA1100154	AC CN L35-L, 276 mm canister height 38 mm diameter	20,00 €
27	AC CN L50-L	H-AAA1100155	AC CN L50-L, 276 mm canister height 38 mm diameter	20,00 €
28	AC CN XL20-L	H-AAA1100156	AC CN XL20-L, 276 mm canister height 63 mm diameter	20,00 €

# AC Medical Family Pricelist




# Auguste Cryogenics

Pos	Model	Part No.	Description	Price
29	AC CN XL25-L	H-AAA1100157	AC CN XL25-L, 276 mm canister height 63 mm diameter	20,00 €
30	AC CN XL30-L	H-AAA1100158	AC CN XL30-L, 276 mm canister height 63 mm diameter	20,00 €
31	AC CN XL35-L	H-AAA1100159	AC CN XL35-L, 276 mm canister height 63 mm diameter	20,00 €
32	AC CN XL50-L	H-AAA1100160	AC CN XL50-L, 276 mm canister height 63 mm diameter	20,00 €
33	AC PG 2L(S1)	H-AAB1100166	AC PG 2L(S1), plug - L (30 mm) - 2 L	20,00 €
34	AC PG 2L(S2)	H-AAB1100167	AC PG 2L(S2), plug - L (35 mm) - 2 L	20,00 €
35	AC PG 3/6/10/13L	H-AAB1100168	AC PG 3/6/10/13L, plug - L -3, 6, 10 & 13 L	21,00 €
36	AC LL L	H-AAB1100178	AC LL L, black lockable lid - L	21,00 €
37	AC PG 15/20/25/30/35/50L	H-AAB1100202	AC PG 15/20/25/30/35/50L, plug - L 15, 20, 25, 30 ,35 & 50 L	21,00 €
38	AC CN 2XL35-S	H-AAA1100147	AC CN 2XL35-S, 120 mm canister height 97 mm diameter	22,00 €
39	AC CN 2XL50-S	H-AAA1100149	AC CN 2XL50-S, 120 mm canister height 97 mm diameter	22,00 €
40	AC LL XL	H-AAB1100179	AC LL XL, black lockable lid - XL	23,00 €
41	AC PG 6/10XL	H-AAB1100169	AC PG 6/10XL, plug - XL - 6 & 10 L	25,00 €
42	AC PG 15/20/25/30/35XL	H-AAB1100170	AC PG 15/20/25/30/35XL, plug - XL - 15, 25, 30, 35, 30 & 35 L	25,00 €
43	AC PG 50XL	H-AAB1100171	AC PG 50XL, plug - XL - 50 L	25,00 €
44	AC CN 2XL30-L	H-AAA1100161	AC CN 2XL30-L, 276 mm canister height 97 mm diameter	30,00 €
45	AC CN 2XL35-L	H-AAA1100162	AC CN 2XL35-L, 276 mm canister height 97 mm diameter	30,00 €
46	AC CN 2XL47-L	H-AAA1100148	AC CN 2XL47-L, Ø97, height 276mm canisters (2 goblets included)	30,00 €
47	AC CN 2XL50-L	H-AAA1100165	AC CN 2XL50-L, 276 mm canister height 97 mm diameter	30,00 €
48	AC PG 10 2XL	H-AAB1100172	AC PG 10 2XL, plug - 2XL - 10 L	34,00 €
49	AC PG 15/20/25/30/35 2XL	H-AAB1100173	AC PG 15/20/25/30/35 2XL, plug - 2XL - 15, 20, 25, 30, 35 L	34,00 €
50	AC PG 47/50 2XL	H-AAB1100174	AC PG 47/50 2XL, plug - 2XL - 47 & 50 L	35,00 €
51	AC LL 2XL	H-AAB1100180	AC LL 2XL, black lockable lid -2XL	35,00 €
52	AC LL 3XL	H-AAB1100181	AC LL 3XL, black lockable lid -3XL	49,00 €
53	AC PG 65/95/115/140/175 3XL	H-AAB1100175	AC PG 65/95/115/140/175 3XL, plug - 3 XL - 175 L	51,00 €
54	AC Cryo gloves 32	H-AAD1100012	AC Cryo gloves 32, 32cm	66,00 €



Pos	Model	Part No.	Description	Price
55	AC Cryo gloves 45	H-AAD1100013	AC Cryo gloves 45, 45cm	72,00 €
56	AC Cryo gloves 72	H-AAD1100014	AC Cryo gloves 72, 72cm	120,00 €
57	AC Multi Port	H-LCB1100388	AC Multi Port, AC multi port connection 6 * G3 / 8	120,00 €
58	AC L-AF	H-LCA1200331	AC L-AF, Flansch für 50mm Halsdurchmesser	160,00 €
59	AC L-AS	H-LCA1200332	AC L-AS, Anschlussstück für Entnahmeheber. 50mm Halsdurchmesser Sicherheitsventil 0,5Bar , Gasentnahme Flüssigkeitsentnahme Manometer	900,00 €
60	AC Roler base small	H-AAD1100016	AC Roler base small , for up to 50 liter vessels	325,00 €
61	AC roler base large	H-AAA1100125	AC roler base large, for 65-175L vessels	360,00 €
62	AC Cryo Low	H-CPB1100387	AC Cryo Low, Low liquid nitrogen level detection	1.500,00 €
63	AC Cryofill	H-CPA1100360	AC Cryofill, multifunctional automatic filling unit for "Pressure differential" or "4 point level"	7.400,00 €
64	AC CN 2XL47-S	H-AAA1100361	AC CN 2XL47-S, Ø72, height 276mm canisters (2 goblets included)	22,00 €
65	AC R25-6	H-AAA1100362	AC R25-6, Rack for 5*5 2ml cryobox; 6 cryoboxes per rack suitable for YDS-50B-125 (include cryobox)	80,00 €
66	AC RR6	H-AAA1100363	AC RR6, Round racks with cryoboxes (hold 37 vials per cryobox) 5 cryoboxes per rack ( cryobox included); suitable for YDS-47-127-6Y	90,00 €
67	AC R25-5	H-AAA1100364	AC R25-5, Rack for 5*5 2ml cryobox; 5 cryoboxes per rack suitable for YDS-35-125, YDS-47-127 (include cryobox)	70,00 €
68	AC CN 2XL30-S	H-AAA1100365	AC CN 2XL30-S, Canisters (Ø97, height 120mm)	30,00 €
69	AC R25-4	H-AAA1100366	AC R25-4, Rack for 5*5 2ml cryobox; 4 cryoboxes per rack suitable for YDS-30-125\YDS-10-125 (include cryobox)	65,00 €
70	AC R81-6	H-AAA1100367	AC R81-6, Rack for 9*9 2ml cryobox; 6 cryoboxes per rack suitable for YDS-115-216; (include cryobox)	110,00 €
71	AC R81-8	H-AAA1100368	AC R81-8, Rack for 9*9 2ml cryobox; 8 cryoboxes per rack suitable for YDS-145-216; (include cryobox)	125,00 €
72	AC R81-1'0	H-AAA1100369	AC R81-1'0, Rack for 9*9 2ml cryobox; 8 cryoboxes per rack suitable for YDS-175-216; (include cryobox)	150,00 €
73	AC R81-4	H-AAA1100370	AC R81-4, Rack for 9*9 2ml cryobox; 4 cryoboxes per rack suitable for YDS-65-216; (include cryobox)	90,00 €
74	AC R81-5	H-AAA1100371	AC R81-5, Rack for 9*9 2ml cryobox; 5 cryoboxes per rack suitable for YDS-95-216, YDS-25H-216; (include cryobox)	100,00 €
75	AC CN DS10-L	H-AAA1100372	AC CN DS10-L, Ø97, height 276mm canisters (2 goblets included) suitable for YDS-10H-125	40,00 €
76	AC CN XL25-S	H-AAA1100373	AC CN XL25-S, Canisters Ø63, height 120mm	20,00 €

# AC Medical Family Pricelist



# Auguste Cryogenics

Pos	Model	Part No.	Description	Price
77	AC CN XL15-S	H-AAA1100374	AC CN XL15-S, Canisters (Ø63, height 120mm)	20,00 €
78	AC CN XL15-L	H-AAA1100375	AC CN XL15-L, Ø63, height 276mm canisters (2 goblets included)	22,00 €
79	AC CN DS3-L	H-AAA1100377	AC CN DS3-L, Canister (Ø38, height 120mm)	22,00 €
80	AC CN 2XL+35-L	H-AAA1100378	AC CN 2XL+35-L, Canisters Ø67, height 276mm	22,00 €
81	AC CN 2XL+47-S	H-AAA1100379	AC CN 2XL+47-S, Ø97, height 120mm canisters	20,00 €
82	AC Cryobox 100	H-AAB1100389	AC Cryobox 100, Hold 100 vials (1.8ml/ 2.0ml; internally threaded)	20,00 €